## Urine Test Strip 5+ Leuko

**Test strips for rapid determination of blood, urobilinogen, protein, nitrite, glucose, and leukocytes in urine**

**Use**

Screening test for detection of diabetes, liver diseases, biliary and hepatic obstructions, hemolytic diseases and diseases of kidney and urinary tract.

**For professional use only.**

**Instructions for use**

Dip the test strip for approximately 1 second into the fresh urine. Draw it across the rim of the container to remove excess urine. After 30 to 60 seconds (leukocyte test field after 60–120 seconds) compare the test strip with the colour scale. The best time for comparison is after 30 seconds. Colour changes that take place after more than 2 minutes are of no significance. When tested the urine should not be older than 2 hours.

**Principle**

**Blood:** The detection is based on the pseudoperoxidative activity of hemoglobin and myoglobin, which catalyse the oxidation of an indicator by an organic hydroperoxide producing a green colour.

**Urobilinogen:** The test paper contains a stable diazonium salt forming a reddish azo compound with urobilinogen.

**Protein:** The test is based on the “protein error” principle of indicators. The test zone is buffered to a constant pH value and changes colour from yellow to greenish blue in the presence of albumin. Other proteins are indicated with less sensitivity.

**Nitrite:** Microorganisms, which are able to reduce nitrate to nitrite, are indicated indirectly by this test. The principle of Griess reagent is the basis of this test. The test paper contains an amine and a coupling component. A red coloured azo compound is formed by diazotisation and subsequent coupling.

**Glucose:** The detection is based on the glucoseoxidase-peroxidase-chromogen reaction. Apart from glucose, no other compound in urine is known to give a positive reaction.

**Leukocytes:** The test is based on the esterase activity of granulocytes. This enzyme splits carboxylic acid ester. The alcohol constituent released reacts with a diazo salt producing a violet colour.

**Evaluation – Sources of Error**

**Blood:** The minimum sensitivity of the test strip is 5 erythrocytes/µL urine corresponding to approx. 0.015 mg hemoglobin/dL urine. Intact erythrocytes are indicated by flecky discolorations of the test field. The colour fields correspond to the following values:

0 (negative), ca. 5–10, ca. 50, ca. 250 Ery/µL resp.

hemoglobin concentration out of ca. 10, ca. 50, ca. 250 Ery/µL

Normal concentrations of ascorbic acid (< 40 mg/dL) do not influence the test results. Falsely positive reactions can be produced by a residue of peroxide containing cleansing agents.

**Urobilinogen:** In dependence upon the urine colour 0.5 to 1 mg urobilinogen/dL urine are indicated. 1 mg/dL is considered to be the normal excretion rate. Higher values are pathological. A complete absence of urobilinogen in the urine, which is likewise pathological, cannot be indicated by the strips. The colour fields correspond to the following urobilinogen concentrations:

norm. (normal), 2, 4, 8, 12 mg/dL or

norm. (normal), 35, 70, 140, 200 µmol/L

The test will be inhibited by higher concentrations of formaldehyde. Exposure of the urine to light for a longer period of time may lead to lowered or falsely negative results. Too high or falsely positive results can be caused by the presence of diagnostic or therapeutic dyes in the urine. Larger amounts of bilirubin produce a yellow colouration.

**Protein:** The minimum sensitivity of the test strip is 10 mg protein/dL urine. The colour fields correspond to the following ranges of albumin concentrations:

negative, 30, 100, 500 mg/dL or

negative, 0.3, 1.0, 5.0 g/L

Falsely positive results are possible in alkaline urine samples (pH > 9), after infusions with polyvinylpyrrolidone (blood substitute), after intake of medicaments containing quinine and also by disinfectant residues in the urine sampling vessel. The protein colouration may be masked by the presence of medical dyes (e.g. methylene blue) or beetroot pigments.

**Nitrite:** The test detects concentrations from 0.05 mg nitrite/dL urine. Every pink colour indicates a bacterial infection of the urinary tract. The colour intensity depends only on the nitrite concentration, but does not provide information about the extent of the infection. A negative result does not preclude an infection of the urinary tract, if bacteria which cannot produce nitrite are present. Falsely negative results can be produced by high doses of ascorbic acid, by antibiotics therapy and by very low nitrate concentrations in urine as the result of low nitrate diet or strong dilution (diuresis). Falsely positive results can be caused by the presence of diagnostic or therapeutic dyes in the urine.

**Glucose:** Pathological glucose concentrations are indicated by a colour change from green to bluish green. Yellow or greenish test fields should be considered negative or normal. The colour fields correspond to the following ranges of glucose concentrations:

neg. (yellow), neg. or normal (greenish), 50, 150, 500, ≥ 1000 mg/dL or

neg. (yellow), neg. or normal (greenish), 2.8, 8.3, 27.8, ≥ 55.5 mmol/L

The influence of ascorbic acid (vitamin C) has been largely eliminated. An inhibitory effect is produced by glutaric acid. Falsely positive reactions can be produced by a residue of peroxide containing cleansing agents.

**Leukocytes:** The test records values starting from approx. 10 leukocytes/µL urine. Changes in colour that can not be assigned to the negative reference field and faint violet colours after 120 seconds must be evaluated as positive. The colour reference fields correspond to the following leukocyte concentrations:

negative (normal), 25, 75, 500 leukocytes/µL

A weakened reaction can be expected in the case of proteinuria of over 500 mg/dL and a glucose concentration of over 2 g/dL as well as in the case of patients taking preparations containing cephalixin and gentamycin. Bacteria, trichomonads and erythrocytes do not react with this test. Formaldehyde (as a preservative) can result in a false positive reaction. Boric acid used as preservative decreases the sensitivity of the reaction. Excretion of bilirubin, nitrofruantoin or other strongly-coloured compounds may disguise the colour of the reaction. Tests with female patients have shown that vaginal discharge can cause a false positive reaction.

**Reactive ingredients**

(quantity resp. activity/cm² at time of impregnation)

<b>Blood:</b>		<b>Nitrite:</b>		<b>Leukocytes:</b>	
tetramethylbenzidine	31 µg	sulfanilic acid	95 µg	carboxylic acid ester	16 µg
cumene hydroperoxide	315 µg	quinoline derivative	37 µg	diazonium salt	14 µg

<b>Urobilinogen:</b>		<b>Glucose:</b>	
diazonium salt	75 µg	glucose oxidase	7 U
		peroxidase	1 U
<b>Protein:</b>		tetramethylbenzidine	96 µg
tetrabromophenol blue	10 µg		

**Quality Control in professional use**

The performance of the test strips should be confirmed by use of positive and negative samples / control solutions. Positive and negative controls should be analyzed once a day, whenever a new bottle of strips is opened, whenever a new lot of strips is started, and every 30 days to check storage conditions. Each laboratory should establish its own goals for adequate standards of performance, and should question handling and testing procedures if these standards are not met.

**Directions**

In any case, in order to establish a final diagnosis and prescribe an appropriate therapy, the results obtained with test strips should be verified with other medical results.

The effect of medicaments or their metabolic products on the test is not known in all cases. In case of doubt it is recommended not to take the medicaments and then repeat the test. Only use well washed and clean vessels for urine collection. The presence of usual urine preservatives will not affect the test results.

Remove only as many test strips as are required, and reseal the container immediately after use. Do not touch the test pads. Avoid exposing the strips to sunlight and moisture. Store the container below + 30 °C in a dry place. The test strips are stable, when stored properly up to the date of expiry indicated.

The caps contain a non-poisonous and harmless desiccant. In case this desiccant is swallowed accidentally, then drink plenty of water.

Explanation of symbols can be found in the package insert.

Disposal: Please dispose all used test strips in accordance with your local laws and regulations.

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**DE Healthcare Ltd**  
**Unit 9, Kingsthorpe Business Centre,**  
**Studland Road, Kingsthorpe,**  
**Northampton NN2 6NE**

Distributed by/Distribuido por/Distribut per/Ventrib distributo/Distribuito del:

And in the USA/Y en EEUU/Et aux États-Unis/Und in USA/E in USA:

**Denver, PA 17517 USA**

**DE Healthcare Products**  
DEDICATED TO EXCELLENCE

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